

## Plate Electrodes Copper Plates

## EM2370-001, EM2371-001 & EM2372-001

## **Description:**

Plates of copper for experiments in electrolytes or for conductivity experiments.

Each plate is fitted with a simple screw terminal for the attachment of a wire.

EM2370-001	Length: 100mm	Width: 50mm	Thickness: 1.5mm
EM2371-001	Length: 100mm	Width: 50mm	Thickness: 3.0mm
EM2372-001	Length: 75mm	Width: 25mm	Thickness: 1.5mm

Designed and manufactured in Australia