

Electrostatic Rubbing Cloths

Cotton flannelette & Wool



EM1772-001

Description:

A set of 2x rubbing cloths, cotton flannelette and wool. 200mm square for rubbing on electrostatic rods for the creation of electrostatic charges.

Specifications:

- 1 x Cotton flannelette cloth, edges overlapped to avoid fraying. Approx 200mm square
- 1 x Wool cloth, edges overlapped to avoid fraying. Approx 200mm square

2 x Cloths	Length: 200mm	Width: 200mm	Weight: 50g
------------	---------------	--------------	-------------

Designed and manufactured in Australia

INDUSTRIAL EQUIPMENT & CONTROL PTY.LTD.

61-65 McClure St. Thornbury. 3071 Melbourne. Australia

Tel: 61 (0)3 9497 2555

Email: iec@iecpl.com.au

www.iecpl.com.au

1