

Electrostatic Rod

Aluminium / Acrylic



EM3550-001

Description:

Special electrostatic rod which is made with one half made from Aluminium metal and the other half made from Acrylic plastic.

Length: 300mm	Diameter: 10mm	Weight: 45g
---------------	----------------	-------------

Designed and Manufactured in Australia

INDUSTRIAL EQUIPMENT & CONTROL PTY.LTD.

61-65 McClure St. Thornbury. 3071 Melbourne. Australia

Tel: 61 (0)3 9497 2555 Email: iec@iecpl.com.au www.iecpl.com.au