

# Weights

## Slotted, Brass, Set of 5



**MF4304-001**

### Description:

The IEC 'Slotted Brass Weights' are solid brass discs with a shaped slot that permits captive engagement on a 5mm diameter rod in the form of a 'carrier' with a strong hook fitted permanently into the top. The set of 5 weights does not include a carrier.

The set is: 1 x 10g, 2 x 20g, 1 x 50g, 1 x 100g to make total weight of 200g.

Each weight is calibrated to an accuracy is better than 1%.

Several types of slotted weights are manufactured by IEC.

Diameter: 30mm	Thickness: Various	Weight: 200g
----------------	--------------------	--------------

Designed and Manufactured in Australia

**INDUSTRIAL EQUIPMENT & CONTROL PTY.LTD.**

61-65 McClure St. Thornbury. 3071 Melbourne. Australia

Tel: 61 (0)3 9497 2555

Email: [iec@iecpl.com.au](mailto:iec@iecpl.com.au)

[www.iecpl.com.au](http://www.iecpl.com.au)

1