

Weights

Slotted, Brass, Set of 3



MF4304-101

Description:

The IEC 'Slotted Brass Weights' are solid brass discs with a shaped slot that permits captive engagement on a 5mm diameter rod in the form of a 'carrier' with a strong hook fitted permanently into the top. The set of 3 weights does not include a carrier.

The set is: 1 x 5g, 1 x 10g, 1 x 20g to make total weight of 35g.

Each weight is calibrated to an accuracy is better than 1%.

Several types of slotted weights are manufactured by IEC.

Diameter: 30mm	Thickness: Various	Weight: 35g
----------------	--------------------	-------------

Designed and Manufactured in Australia