

# Weights, Slotted

## Set of 9, 100g Weights on Carrier



**MF4301-001**

### Description:

The IEC 'Slotted Brass Weights' are solid brass discs with a shaped slot that permits captive engagement on a 5mm diameter rod in the form of a 'carrier' with a strong hook fitted permanently into the top. The set of 9 weights fitted to a 100g carrier makes a total weight of 1000g.

Each weight in the set is marked 100g together with the brass carrier with hook.

Each weight and the carrier is calibrated to an accuracy is better than 1%.

Several types of slotted weights are manufactured by IEC.

Diameter: 30mm	Height: 217mm	Weight: 1.0kg
----------------	---------------	---------------

Designed and manufactured in Australia