

Weights, Slotted

Set of 9, 50g Weights on Carrier



MF4300-001

Description:

The IEC '**Slotted Brass Weights**' are solid brass discs with a shaped slot that permits captive engagement on a 5mm diameter rod in the form of a 'carrier' with a strong hook fitted permanently into the top. The set of 9 weights fitted to a 50g carrier makes a total weight of 500g. Each weight in the set is marked 50g together with the brass carrier with hook.

Each weight and the carrier is calibrated to an accuracy is better than 1%.

Several types of slotted weights are manufactured by IEC.

Diameter: 30mm	Height: 125mm	Weight: 500g
----------------	---------------	--------------

Designed and manufactured in Australia

INDUSTRIAL EQUIPMENT & CONTROL PTY.LTD.

61-65 McClure St. Thornbury. 3071 Melbourne. Australia

Tel: 61 (0)3 9497 2555

Email: iec@iecipl.com.au

www.iecipl.com.au

1