

Light Source

Very High Power



HL2016-001 100W Halogen Lamp

Description:

The IEC High Power Light Source is used as a convenient source of very bright white light.

It can be viewed from either of two directions and each viewing hole has a swinging shutter.

The light source can be supported by clamping to the rod protruding from the rear of the body.

Power connection to the lamp is by 4mm banana sockets at the lower edge of the housing.

The lamp runs from 12V.AC or DC and draws about 8.3 amps

The lamp is Quartz Halogen, 2 pin, 12V. 100W Supplied complete with lamp, cable and stackable banana plugs.

Lamp Type: 12v 100W 2 pin, quartz halogen, lateral filament, PA2043-001.

Length: 76mm Width: 76mm Height: 153mm Rod: 150mm Weight: 5.6kg

Designed and manufactured in Australia

INDUSTRIAL EQUIPMENT & CONTROL PTY.LTD.

27-Oct-21