

Calorimeter Box & Lid

Styrene Foam Insulation



HL0830-001 Insulation Box & Lid

Description:

The 'IEC' **calorimeter insulating box** and lid is designed to accept the standard 50mm diameter x 75 deep calorimeter cup. It is moulded from highly insulating polystyrene foam and is attractive and very durable. The lid engages over the mouth of the lower section and firmly fits in place and is not easily dislodged.

The lid has a central hole for thermometer, a hole to one side for a stirrer to pass into the calorimeter and holes can easily be created to pass over the 4mm screw terminals when used with a Joule's Calorimeter.

Height: 120mm	Diameter: 72mm	Weight: 50g
---------------	----------------	-------------

Designed and manufactured in Australia