

Magnetic Needle on Stand

75mm Long Needle



EM2070-001

Description:

The magnetic needle is fitted with a brass socket which rests on the vertical point of the stand so that it is very free to rotate in the horizontal plane. One end of the needle is marked to indicate the end that points to the magnetic north. This end is actually the south pole of the needle. When a weak magnetic field is placed close to the needle, it will deflect from rest.

Diameter: 75mm	Height: 150mm	Weight: 80g
----------------	---------------	-------------

Designed and manufactured in Australia