

Scale Pan

Nylon, with Hook



MF3723-001

Description:

The IEC 'Scale Pan with Hook' is a low cost pan for use in school laboratories for many purposes. It is 100mm diameter and 12mm deep strong plastic moulding with integral hinged side flaps that converge to form a hook and supports.

The picture shows the pan with the side flaps converged to form the hook.

The pan can be stored flat. The arms are pivoted and brought together to form the pan hangers.

Diameter: 100mm	Height: 12mm	Weight: 10g
-----------------	--------------	-------------

Designed and Manufactured in Australia

INDUSTRIAL EQUIPMENT & CONTROL PTY.LTD.

61-65 McClure St. Thornbury. 3071 Melbourne. Australia

Tel: 61 (0)3 9497 2555

Email: iec@iecpl.com.au

www.iecpl.com.au

1