

# Scale Pan

## Metal with Hook



**MF3721-001 Brass**



**MF3722-001 Aluminium**

### Description:

The IEC ‘**Scale Pan with Hook**’ is the classic scale pan used in school laboratories for many purposes. It is 75mm diameter and 12mm deep with 3 holes around the edge to attach 3 cords or chains. These cords or chains connect to a single hook which is then used to support the pan. They are used in the laboratory for general purpose weighing or for carrying weights for experiments, for example, Hooke’s Law. These pans are available in brass or aluminium.

IEC manufactures also a fully moulded nylon scale pan MF3723-001 which is very useful and easy to use and store.

|                |              |                                      |
|----------------|--------------|--------------------------------------|
| Diameter: 75mm | Length: 12mm | Weight: 10g (Aluminium), 30g (Brass) |
|----------------|--------------|--------------------------------------|

Designed and Manufactured in Australia