

Switch

SPDT / SPST, Toggle

**LB3960-050 SPDT Switch****LB3960-060 SPST Switch**

Description:

The IEC Toggle Switches are compact and strong units that can be used in general electrical experiments. Sometimes they are standard components from the various excellent and robust IEC Electricity Teaching Kit. e.g. EM1763-001 & EM1764-01

This sheet covers 2 types:

- SPST (Single Pole Single Throw). This means that only one circuit can be controlled and that the control is only ON and OFF. This is the simplest form of switch.
- DPST (Double Pole Single Throw) This means that 2 circuits are controlled, but the switching in each is only ON and OFF for both circuits at the same time.

The standard symbol for each switch is printed clearly in the housing below the switch.

For export to Europe, we can supply these and many other electrical items with 4mm safety banana sockets.

Note:

Many other parts of the Electricity Kit EM1763-001 and the Electronics Kit EM1764-001 are useful in the classroom as individual components. Please refer to the PA (Parts & Accessories) section of our listing for parts from these kits.

Length: 75mm	Width: 35mm	Height: 35mm	Weight: 100g
--------------	-------------	--------------	--------------

Designed and manufactured in Australia

INDUSTRIAL EQUIPMENT & CONTROL PTY.LTD.

61-65 McClure St. Thornbury. 3071 Melbourne. Australia

Tel: 61 (0)3 9497 2555 Fax: 61 (0)3 9497 2166 www.iecpl.com.au