


Plate Electrodes

Iron Plates



EM2380-001 & EM2381-001

Description:

Plates of iron for experiments in electrolytes or for conductivity experiments.

Each plate is fitted with a simple screw terminal for the attachment of a wire.

EM2380-001	Length: 100mm	Width: 50mm	Thickness: 1.5mm
EM2381-001	Length: 75mm	Width: 25mm	Thickness: 1.5mm

Designed and manufactured in Australia

INDUSTRIAL EQUIPMENT & CONTROL PTY.LTD.

61-65 McClure St. Thornbury, 3071 Melbourne, Australia

Tel: 61 (0)3 9497 2555 Email: iec@iecpl.com.au www.iecpl.com.au

1