


# Plate Electrodes

## Lead Plates



### EM2390-001, EM2391-001 & EM2392-001

#### Description:

Plates of lead for experiments in electrolytes or for conductivity experiments.

Each plate is fitted with a simple screw terminal for the attachment of a wire.

EM2390-001	Length: 100mm	Width: 50mm	Thickness: 1.5mm
EM2391-001	Length: 100mm	Width: 50mm	Thickness: 3.0mm
EM2392-001	Length: 75mm	Width: 25mm	Thickness: 1.5mm

Designed and manufactured in Australia

**INDUSTRIAL EQUIPMENT & CONTROL PTY.LTD.**

61-65 McClure St. Thornbury, 3071 Melbourne, Australia

Tel: 61 (0)3 9497 2555 Email: [iec@iecpl.com.au](mailto:iec@iecpl.com.au) [www.iecpl.com.au](http://www.iecpl.com.au)

1