



Portable Biological Microscope

4x, 10x, 40x, LED Illumination

- **4x, 10x, 40x Objective Lens**
- **WF 10x/20mm Eyepieces**
- **Upper & Lower LED Light**
- **3 AAA Batteries**
- **Optional DC Power Plug**

The MS1012 is a versatile and compact microscope designed for easy use in various settings. It features three objective lenses with magnifications of 4x, 10x, and 40x, allowing users to observe specimens at different levels of detail.

Applications

The MS1012 Portable Biological Microscope offers flexibility, convenience, and versatility, making it suitable for students, hobbyists, or professionals who require a portable microscope for biological observations and research.



Product Specifications

Head	Monocular
Eyepiece	Wide Field 10x/20mm
Magnification	4x, 10x, 40x Objective Lens
Illumination	Upper and Lower LED light source
Batteries	3 x AA Batteries
Power Supply	Optional DC Plug Pack Available

