


Proof Plane

Flat Disc or Sphere



EM2520-001 & EM2521-001

Description:

The classic Proof Plane is a 20mm diameter metal disc or metal sphere attached to a insulating handle. It is used to collect electrostatic charges from devices and to deposit them on other devices.

Design can be either brass or aluminium alloy and the handle can be either acrylic or polystyrene. Performance is exactly the same.

EM2520-001	Length: 140mm	Diameter: 20mm	Weight:: 22g
EM2521-001	Length: 140mm	Diameter: 20mm	Weight:: 40g

Designed and manufactured in Australia